




The Audicus Mini 2 - Rechargeable

Fitting Range



Technical Data

	2cc Coupler
OSPL90 Maximum	114 dB SPL
HF Average OSPL90	108 dB SPL
Peak Gain	50 dB
HFA Full-On-Gain	45 dB
Reference Test Gain	31 dB
Frequency Range	100 - 9100 Hz
THD at 500 Hz	<2%
800 Hz	<3%
1600 Hz	<2%
Equivalent Input Noise	16 dB SPL
Latency	<15 ms
Battery runtime	Up to 28 h

2cc Coupler

